




# C-Slide

*The New Name  
for Countess Slides*




For use with Automated Cell Counters



# C-Slide

## Cell Counter Chamber Slides

- No need for coverslips
- Disposable, eliminates the need for washing
- Reduce exposure to infectious samples and risk of contamination
- Guaranteed reproducibility and reliability
- Made of sturdy quartz grade optical plastic



NanoEnTek

The developer & manufacturer  
of Countess and Countess slide.

## Specifications

Sample volume : 10  $\mu$ L

Staining method : Trypan blue stain

Dimension : 75 (L) x 25 (W) x 1.8 (H) mm

## For use with Automated Cell Counters



Countess



Countess II



## Ordering information

Distributed in the UK & Ireland by:



Visit: [www.bioscience.co.uk](http://www.bioscience.co.uk)

Call: 01223 316855

Email: [sales@bioscience.co.uk](mailto:sales@bioscience.co.uk)

Catalog Number	Product Name
CS-050	C-Slide, Cell Counter Chamber Slides 50 slides (100 tests)

**NanoEnTek**

e-mail [sales@nanoentek.com](mailto:sales@nanoentek.com)  
website [www.nanoentek.com](http://www.nanoentek.com)

**NanoEnTek, Inc.**

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,  
Gyeonggi-do, 18531, Korea  
Tel: +82-2-6220-7940 / Fax: +82-2-6220-7999

**NanoEnTek America, Inc.**

240 Bear Hill Road, Suite 101, Waltham,  
MA 02451, USA  
Tel: +1-781-472-2558 / Fax: +1-781-790-5649

**EC REP**

MT Promedt Consulting GmbH  
Altenhofstrasse 80,  
66386 St. Ingbert, Germany